

Transport Equipment

Single-use transport equipment is widely used in biopharmaceutical processes for the storage and transportation of liquids. LePure provides complete system transportation and storage solutions to support the needs of liquid handling in the production process, meeting transportation and handling. EZ-TRANS series transfer tools, specially designed for 2D and 3D reservoirs, suitable for different volumes. Different shapes of reservoir bags, equipped with different transfer tools, can achieve flexible transfer.



Typical applications

- Storage
- Transfer
- Mix

Advantage

- Small footprint
- Foldable and stackable
- Easy trans
- Economically applicable
- Diversified capacity
-

Technical specification

Foldable tank

Article No	COMB0250E	IBC500	COMB01000E
Material	PP	HDPE	PP
External dimensions(with lid, mm)	809x609x833	1210x810x932	1155x1155x1132
External dimensions(without lid, mm)	800x600x810	1200x800x910	1145x1145x1110
Internal dimensions(mm)	740x540x660	1110x710x690	1065x1045x940
Folded height(mm)	403	430	467
Floor space(m ²)	0.5	0.96	1.33
Container capacity(kg)	350	650	1400
Weight(kg)	20	67	89
Max load stack, dynamic(Layer)	2	2	2
Max container capacity(L)	264L	520L	1060L
Temperature range	-20°C to +40°C	-20°C to +40°C	-20°C to +40°C

Supporting base	COMB0250E	IBC500	COMB01000E
Supporting base(Material)	SS304	SS304	SS304
Supporting base(Dimensions mm)	830x630x230	1250x850x230	1190x1190x230
Supporting base(Weight kg)	12.5	38	42.5

2D bag plastic tray

Article No	LTR005	LTR010	LTR020	LTR030	LTR050
Material	PP	PP	PP	PP	PP
Size	530*330*125	590*440*145	780*480*165	780*530*185	930*630*200
Supporting base	DOLLYCL005	DOLLYCL010	DOLLYCL020	DOLLYCL030	DOLLYCL050
Supporting base size	550*360	610*480	800*500	800*550	950*680
Brand	LePure	LePure	LePure	LePure	LePure
Delivery	2-4 weeks	2-4 weeks	2-4 weeks	2-4 weeks	2-4 weeks

Stainless steel barrel

Article No	BIN028	BIN050	BIN100	BIN200	BIN303	BIN378
Material	304	304	304	304	304	304
Size	R303*H500	R353*H694	R454*H760	R554*H910	R606*H1180	R706*H1120
Supporting base	DOLLY028	DOLLY057	DOLLY100	DOLLY200	DOLLY303	DOLLY378
Supporting base size	R350	R380	R490	R580	R630	R766
Brand	LePure	LePure	LePure	LePure	LePure	LePure
Delivery	2-4 weeks	2-4 weeks	2-4 weeks	2-4 weeks	2-4 weeks	2-4 weeks

Round plastic bucket

Article No	Size	Supporting base	Supporting base size	Material	Brand	Delivery
(C)LBK019-1	R279*H384	DOLLY019	R320	304	LePure	2-4 weeks
(C)LBK028-1	R316*H457	DOLLY028	R350			
(C)LBK038-1	R340*H518	DOLLY038	R380			
(C)LBK057-1	R346*H702	DOLLY057	R490			
(C)LBK0113-1	R467*H759	DOLLY113	R580			
(C)LBK0208-1	R559*H918	DOLLY208	R630			
(C)LBK0303-1	R608*H1221	DOLLY303	R630			
(C)LBK0378-1	R721*H1124	DOLLY378	R766			
(C)LBK0568-1	R803*H1251	DOLLY568	R840			
(C)LBK0757-1	R932*H1295	DOLLY757	R940			

Article No	SV50517.04	SV50517.05	SV50517.06	SV50517.06	SV50517.09	SV50517.10
Material	LLDPE	LLDPE	LLDPE	LLDPE	LLDPE	LLDPE
Size	R600*580	R600*760	R600*1140	R600*1140	R600*760	R600*1140
Technical Information	Conical drum, 50L, 1-port 4.5cm (1.75 in.), with fixture	Conical drum, 100L, 1-port 4.5cm (1.75 in.), with fixture	Conical drum, 200L, 1-port 4.5cm (1.75 in.), with fixture	Conical drum, 50L, 1 port 10.2 cm (4 in.), with fixture	Conical drum, 100L, 1 port 10.2 cm (4 in.), with fixture	Conical drum, 200L, 1 port 10.2 cm (4 in.), with fixture

Ordering information

Foldable tank

Name	Article No	Packaging
250L Foldable tank	COMB0250E	NA
500L Foldable tank	IBC500	NA
1000L Foldable tank	COMBO1000E	NA

2D bag plastic tray

Name	Article No	Packaging
5L 2D bag plastic tray	LTR005	1 unit/box
10L 2D bag plastic tray	LTR010	1 unit/box
20L 2D bag plastic tray	LTR020	1 unit/box
30L 2D bag plastic tray	LTR030	1 unit/box
50L 2D bag plastic tray	LTR050	1 unit/box

Stainless steel barrel

Article No	Article No	Packaging
28L Stainless steel barrel	BIN028	NA
50L Stainless steel barrel	BIN050	NA
100L Stainless steel barrel	BIN100	NA
200L Stainless steel barrel	BIN200	NA
303L Stainless steel barrel	BIN303	NA
378L Stainless steel barrel	BIN378	NA

Round plastic bucket

Name	Article No	Packaging
19L Round plastic bucket	(C) LBK019-1	1 unit/box
28L Round plastic bucket	(C) LBK028-1	1 unit/box
38L Round plastic bucket	(C) LBK038-1	1 unit/box
57L Round plastic bucket	(C) LBK057-1	1 unit/box
113L Round plastic bucket	(C) LBK0113-1	1 unit/box
208L Round plastic bucket	(C) LBK0208-1	1 unit/box
303L Round plastic bucket	(C) LBK0303-1	1 unit/box
378L Round plastic bucket	(C) LBK0378-1	1 unit/box
568L Round plastic bucket	(C) LBK0568-1	1 unit/box
757L Round plastic bucket	(C) LBK0757-1	1 unit/box
50L Round plastic bucket	SV50517.04	1 unit/box
100L Round plastic bucket	SV50517.05	1 unit/box
200L Round plastic bucket	SV50517.06	1 unit/box
50L Round plastic bucket	SV50517.08	1 unit/box
100L Round plastic bucket	SV50517.09	1 unit/box
200L Round plastic bucket	SV50517.10	1 unit/box